

How does **MEDIHONEY® / Active *Leptospermum* Honey** help to support the healing of routine and challenging wounds?

Gel vs *Paste Dressings*

- Helps to maintain a moist environment for optimal healing¹⁻⁶
- High osmotic potential pulls fluid from deeper tissues to the wound surface, aiding in autolytic debridement¹⁻⁶
- Low pH helps to lower the pH within the wound environment^{7,8}, which has been shown to have wound healing benefits.⁹
- Is easy-to-use, non-toxic, natural and safe with the potential for extended wear times

Gel



Paste



- 80% Active *Leptospermum* Honey and 20% Natural Gelling agents
- More viscous formulation with increased stability at the site of the wound
- Demonstrated positive clinical outcomes^{1, 5, 6}
- Optimal for superficial to full thickness wounds

- 100% Active *Leptospermum* Honey
- Demonstrated positive clinical outcomes²⁻⁴
- Optimal for deep tunneled or undermined wounds

MEDIHONEY® Ordering Information

Order Code	Description	Packaging unit/Case	HCPCS*	
Gel				
31805	0.5 oz tube	10/box, 4 boxes/case	A4649	
31815	1.5 oz tube	1/box, 12 boxes/case	A4649	
HCS				
Non-adhesive				
31622	2.4" x 2.4"	10/box, 5 boxes/case	A4649	
31644	4.33" x 4.33"	10/box, 5 boxes/case	A4649	
31688	8" x 8"	5/box, 4 boxes/case	N/A	
31612	8" x 12"	2/box, 5 boxes/case	N/A	
Fenestrated - Non-adhesive				
31618	1.8" x 1.8"	10/box, 5 boxes/case	A4649	
Adhesive				
31722	2.8" x 2.8" (4.3" x 4.3" with adhesive border)	10/box, 5 boxes/case	A4649	
31744	4 1/2" x 4 1/2" (6" x 6" with adhesive border)	10/box, 5 boxes/case	A4649	
Calcium Alginate				
31012	3/4" x 12"	5/box, 4 boxes/case	A4649	
31022	2" x 2"	10/box, 10 boxes/case	A4649	
31045	4" x 5"	10/box, 5 boxes/case	A4649	
Honeycolloid™				
Non-adhesive				
31222	2" x 2"	10/box, 10 boxes/case	A4649	
31245	4" x 5"	10/box, 5 boxes/case	A4649	
Adhesive				
31422	2" x 2" (3 1/2" x 3 1/2" with adhesive border)	10/box, 10 boxes/case	A4649	
31445	4 1/2" x 4 1/2" (6" x 6" with adhesive border)	10/box, 5 boxes/case	A4649	
Paste				
31505	0.5 oz tube	10/box, 4 boxes/case	A4649	
31515	1.5 oz tube	1/box, 12 boxes/case	A4649	
31535	3.5 oz tube	1/box, 12 boxes/case	A4649	

*Refer to www.dmeptac.com for the most current HCPCS coding of MEDIHONEY® surgical dressings.

Pair MEDIHONEY® with our super absorbent cover dressing, XTRASORB®. It's osmotic gradient pulls exudate to the back of the dressing and converts it into a gel, locking it away - even under compression!

